Dräger

cs

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Dräger DrugCheck[®] 1200 UTK

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11 Odkazy na literaturu

1. Urine Testing for Drugs of Abuse, National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.
 2. R. C. Baselt, Disposition of Toxic Drugs and Chemicals in Man, 2nd Ed., Biomedical Publications, Davis Ca., 1982.
 3. Federal Register, Department of Health and Human Services Mandatory Guidelines for Workplace Drug Testing Program, 59, 110, 22918-29931 (1994).
 4. Stability of drugs of abuse in urine samples stored at -20 °C. S. Dugan, et. Al. J. Anal. Tox. 18 (7) 391-396 (1994).
 5. Long-term stability of abused drugs and anti-abuse chemotherapeutic agents stored at -20 °C., D.E. Moody, et.al., J. Anal. Tox. 23 (6) 535-540 (1999).
 6. CDC (1987) Universal Precautions for Prevention of Transmission of HIV and Other Bloodborne Infections. MMWR 1988,37:377-388
- Testovací kazeta Dräger DrugCheck 1200 UTK je vyvíjena a vyráběna firmou American Bio Medica Corporation.

Methamphetamines	concentration [ng/mL]
Procaine (Novocaine)	60,000
Trimethobenzamide	20,000
+/- methamphetamine	1000
+ methamphetamine	500
Ranitidine (Zantac)	50,000
(+/-) 3,4-Methylenedioxymethamphetamine (MDMA)	1000

Opiates 300 ng/mL	concentration [ng/mL]
6-Acetylmorphine	500
Codeine	100
Eserine (Physostigmine)	15,000
Ethylmorphine	100
Heroin (Diacetylmorphine)	500
Hydromorphone	2000
Hydrocodone	1250
Morphine	300
Morphine-3-b-D-Glucuronide	75
Nalorphine	500
Norcodeine	35,000
Oxycodone	50,000
Thebaine (Paramorphine)	13,000

Opiates 2000 ng/mL	concentration [ng/mL]
6-Acetylmorphine	1000
Codeine	800
Ethylmorphine	400
Heroin (Diacetylmorphine)	10,000
Hydromorphone	2000
Hydrocodone	5000
Morphine	1600
Morphine-3-b-D-Glucuronide	2000
Oxycodone	75,000
Thebaine (Paramorphine)	26,00

Oxycodone	concentration [ng/mL]
6-Acetyl codeine	25,000
6-acetylmorphine	75,000
Codeine	12,500
Dihydrocodeine	3125
Hydromorphone	2500
Hydrocodone	625
Morphine	6250
Noroxycodone	50,000
Oxycodone	100
Oxymorphone	100
Thebaine	25,000

Phencyclidine (PCP)	concentration [ng/mL]
Phencyclidine	25
4-Hydroxy phencyclidine	90

Phencyclidine (PCP)	concentration [ng/mL]
Phencyclidine Morpholine	625

Dräger DrugCheck 1200 UTK PCP also detect high concentrations of the cough suppressant, dextromethorphan. In young children, dextromethorphan overdoses may produce a preliminary positive result for PCP. However, adults ingesting therapeutic dosages of dextromethorphan should not produce a preliminary positive result.

Propoxyphene	concentration [ng/mL]
Propoxyphene	300
Norpropoxyphene	300

THC/Cannabinoids (Tetrahydrocannabinol)	concentration [ng/mL]
Cannabinol	
Efavirenz*	25,000
11-Hydroxy-D9-Tetrahydrocannabinol	5000
11-Nor-D8-Tetrahydrocannabinol-9 Carboxylic Acid	50
11-Nor-D9-Tetrahydrocannabinol-9 Carboxylic Acid	50
11-Nor-D9-Tetrahydrocannabinol-9 Carboxylic Acid Glucuronide	2500
D8-Tetrahydrocannabinol	20,000
D9-Tetrahydrocannabinol	20,000

* Efavirenz is the generic drug found in some HIV treatment medications. Research sources have indicated that it is highly possible false positive results for THC may be observed in patients taking medications which may include Efavirenz.

Tricyclic Antidepressants	concentration [ng/mL]
Amitriptyline	1000
Clomipramine	75,000
Cyclobenzaprine	8000
Cyproheptadine	50,000
Desipramine	1000
Doxepin	5000
Imipramine	1000
Norclomipramine	2500
Nordoxepin	500
Nortriptyline	1000
Promazine	12,500
Protriptyline	2000
Trimipramine	3000

T2 - Sensitivity

Drug	average concentration [ng/mL]
Amphetamines	1000
Barbiturates	300
Benzodiazepines	300
Buprenorphine	12.5
Cocaine	150
Cocaine	300
MDMA	1000
Methadone	300
Methamphetamines	1000
Opiates	300
Opiates	2000
Oxycodone	100
Phencyclidine	25
Propoxyphene	300
THC/Cannabinoids	50
Tricyclic Antidepressants	1000

T3 - Summary

Drug	concentration [ng/mL]	results # pos./10
Amphetamines	500	0/10
	750	2/10
	1000	10/10
	1250	10/10
Barbiturates	150	0/10
	225	2/10
	300	10/10
	375	10/10
Benzodiazepines	150	0/10
	225	2/10
	300	10/10
	375	10/10
Buprenorphine	5	0/10
	10	2/10
	12.5	10/10
	15	10/10
Cocaine 150 ng/mL	75	0/10
	113	2/10
	150	10/10
	187	10/10
Cocaine 300 ng/mL	150	0/10
	225	3/10
	300	10/10
	375	10/10
MDMA	500	0/10
	750	2/10
	1000	10/10
	1250	10/10
Methadone	150	1/10
	225	3/10
	300	10/10
	375	10/10
Methamphetamines	500	0/10
	750	3/10
	1000	10/10
	1250	10/10
Opiates 300 ng/mL	150	0/10
	225	2/10
	300	10/10
	375	10/10
Opiates 2000 ng/mL	1000	0/10
	1250	3/10
	2000	10/10
	2500	10/10
Oxycodone	50	0/10
	75	3/10
	100	10/10
	125	10/10
Phencyclidine	13	0/10
	19	3/10
	25	10/10
	37	10/10
Propoxyphene	150	0/10
	225	3/10
	300	10/10
	375	10/10
THC/Cannabinoids	25	0/10
	38	3/10
	50	10/10
	75	10/10
Tricyclic Antidepressants	500	0/10
	750	2/10
	1000	10/10
	1250	10/10

T4 - Accuracy

drug name	result [ng/mL]	Negative less than 50% of the Cut-off concentration	Near Cut-off negative (between 50% below the Cut-off and the Cut-off concentration)	Near Cut-off positive (between the Cut-off and 50 % above the Cut-off concentration)	High positive (greater than 50% above the Cut-off concentration)	percent agreement
Amphetamines	Positive	0	7	21	19	100 %
	Negative	13	21	0	0	83 %
Barbiturates	Positive	0	3	11	34	100 %
	Negative	16	12	0	0	90 %
Benzodiazepines	Positive	0	7	13	27	100 %
	Negative	15	21	0	0	84 %
Buprenorphine*	Positive	0	6	11	30	100 %
	Negative	13	20	0	0	85 %
Cocaine 150 ng/mL	Positive	0	7	11	33	100 %
	Negative	17	15	0	0	82 %
Cocaine 300 ng/mL	Positive	0	4	15	30	100 %
	Negative	15	17	0	0	89 %
MDMA	Positive	0	9	8	33	100 %
	Negative	6	29	0	0	80 %
Methadone	Positive	0	7	10	30	100 %
	Negative	22	14	0	0	80 %
Methamphetamines	Positive	0	8	14	26	100 %
	Negative	10	26	0	0	82 %
Opiates 300 ng/mL	Positive	0	1	10	50	100 %
	Negative	11	8	0	0	95 %
Opiates 2000 ng/mL	Positive	0	11	12	28	100 %
	Negative	29	19	0	0	81 %
Oxycodone	Positive	0	1	9	37	100 %
	Negative	10	13	0	0	81 %
Phencyclidine	Positive	0	4	7	33	100 %
	Negative	8	14	0	0	72 %
Propoxyphene	Positive	0	13	18	22	100 %
	Negative	4	30	0	0	72 %
THC/Cannabinoids	Positive	0	4	17	27	100 %
	Negative	15	17	0	0	89 %
Tricyclic Antidepressants**	Positive	0	5	15	25	100 %
	Negative	15	17	0	0	86 %

* Confirmation was done with LC/MS

** Confirmation was done with HPLC

T5 - Reproducibility

drug	concentration [ng/mL]	#	results	precision
Amphetamines	500	40	40 neg	>99 %
	750	40	34 neg	>85 %
	1000	40	40 pos	>99 %
	1250	40	40 pos	>99 %
Barbiturates	150	40	40 neg	>99 %
	225	40	36 neg	>90 %
	300	40	40 pos	>99 %
	375	40	40 pos	>99 %
Benzodiazepines	150	40	40 neg	>99%
	225	40	34 neg	>85%
	300	40	40 pos	>99%
	375	40	40 pos	>99%
Buprenorphine	6.3	40	39 neg	>97%
	9.4	40	33 neg	>82%
	12.5	40	40 pos	>99%
	15.6	40	40 pos	>99%
Cocaine 150 ng/mL	75	40	40 neg	>99%
	113	40	33 neg	>82%
	150	40	40 pos	>99%
	187	40	40 pos	>99%
Cocaine 300 ng/mL	150	40	40 neg	>99%
	225	40	34 neg	>85%
	300	40	40 pos	>99%
	375	40	40 pos	>99%
MDMA	500	80	80 neg	>99%
	750	80	65 neg	>81%
	1000	80	80 pos	>99%
	1250	80	80 pos	>99%
Methadone	150	40	39 neg	>97%
	225	40	36 neg	>90%
	300	40	40 pos	>99%
	375	40	40 neg	>99%
Methampheta- mines	500	40	40 neg	>99%
	750	40	36 neg	>90%
	1000	40	40 pos	>99%
	1250	40	40 pos	>99%
Opiates 300 ng/mL	150	40	40 neg	>99%
	225	40	32 neg	>80%
	300	40	40 pos	>99%
	375	40	40 pos	>99%
Opiates 2000 ng/mL	1000	40	40 neg	>99%
	1250	40	31 neg	>77%
	2000	40	40 pos	>99%
	2500	40	40 pos	>99%
Oxycodone	50	80	72 neg	>90%
	75	80	65 pos	>81%
	100	80	80 pos	>99%
	125	80	80 pos	>99%
Phencyclidine	13	40	40 neg	>99%
	19	40	35 neg	>87%
	25	40	40 pos	>99%
	37	40	40 pos	>99%
Propoxyphene	150	40	40 neg	>99%
	225	40	34 neg	>85%
	300	40	40 pos	>99%
	375	40	40 pos	>99%
THC/ Cannabinoids	25	40	40 neg	>99%
	38	40	32 neg	>80%
	50	40	40 pos	>99%
	75	40	40 pos	>99%
Tricyclic Antidepressants	500	40	40 neg	>99%
	750	40	34 neg	>85%
	1000	40	40 pos	>99%
	1250	40	40 pos	>99%

T6 - Cross reactivity

Acebutolol Acetaldehyde Acetaminophen (4-Acetamido-phenol; N-Acetyl-paminophenol) Acetazolamide Acetone 3-(α-acetylbenzyl)-4-hydroxy-coumarin (Warfarin) Acetophenetidin Acetopromazine N-Acetyl-L-cysteine 6-Acetylmorphine (except OPI & OXY) N-Acetylprocainamide (Acedainide) Acetylsalicylic Acid (Aspirin; 2-Acetoxybenzoic Acid) Albumin, standard Albuterol Allobarbitol (5,5-Diallylbarbituric Acid) (except BAR) Allopurinol (4-Hydroxypyrazole (3,4-) Pyrimidine) Alpha-hydroxyriazolam* Alprazolam (except BZO) Alprenolol Amantadine (Adamantan-1-amine) Amcinonide (+) Amethopterin (4-Amino-10-methylfolic acid; Methotrexate; Methylaminopterin) Amikacin Amiloride p-Aminobenzoic Acid 7-Aminoclonazepam 7-Aminoflunitrazepam DL-Aminoglutethimide 7-Aminonitrazepam Amiodarone Amitriptyline (except TCA) Ammonium Chloride Amobarbital (amytal; 5-Ethyl-5-Isoamyl barbituric Acid) (except BAR) Amoxetone Amoxicillin Amphoterin B D-Amphetamine (except AMP) DL-Amphetamine (except AMP) L-Amphetamine (except AMP) Ampicillin D-Amygdalin Aniline Antipyrine (Phenazone) Apomorphine Aprobarbital (except BAR) Aripiprazole Aripiprazole (-) Arterenol [(-)Norepinephrine] L-Ascorbic Acid ASP-PHE-Methyl-Ester (Aspartame) D-Aspartic Acid DL-Aspartic Acid L-Aspartic Acid Astemizole Atenolol Atomoxetine Atropine (Tropinotropate) Atrovastin Azathioprine Baclofen Barbitol (Barbitone; 5,5-Diethylbarbituric acid; Veronal) (except BAR) Barbituric Acid (2,4,6-Trihydroxypyrimidine; Malonylurea) Beclomethasone Beclomethasone Dipropionate Bendroflumethiazide Benzidine (4,4 Diaminobiphenyl) Benicar	Benzylic Acid β-diethylaminoethyl ester Benzocaine (Ethyl-p-Aminobenzoate) Benzoic Acid Benzonatate Benzoylecgonine (except COC) Benzphetamine (α-dimethylphenethylamine) Benzthiazide Bzotropine Methane sulfonate (Benztropine Mesylate) Benzyl alcohol Benzylamine Benzylpiperazine Berberine Betamethasone Bilirubin Bisacodyl Bromazepam (except BZO) 2-Bromo-α-ergocryptine (Bromocriptine mesylate) (+) Brompheniramine (Dexbrompheniramine) (+/-) Brompheniramine Bumetanide Bupivacaine Buprenorphine (except BUP) Bupropion HCL Buspirone Butabarbital (except BAR) Butalbital (except BAR) Butethal (except BAR) Butacaine 2-Butynoic Acid Ethyl Ester (Ethyl-2-Butynoate) Butyrophene Caffeine (1,3,7-Trimethylxanthine) (+/-) Camphor Cannabidiol Cannabinol (except THC) Canrenoic Acid Captopril Carbamazepine Carbamyl-β-methylcholine-chloride (Bethanechol Chloride) Carboplatin (S)-(-)-Carbidopa Carisoprodol Carvedilol Cefaclor Cefadroxil Cefotaxime Cefoxitin Ceftriaxone Cefuroxime Cephalixin Cephaloridine Cephadrine (Cefradine) Cetirizine α-Chloralose Chloramphenicol (Chloromycetin) Chlorcyclizine Chlordiazepoxide (except BZO) 2-(p-Chlorophenoxy)-2-Methylpropionic Acid Ethyl Ester (Clobifate) Chloroquine Chlorothiazide Chlorotrianisene (+)Chlorpheniramine (+/-)Chlorpheniramine Chlorpromazine Chlorpropamide Chlorprothixene Chlorthalidone Chlorzoxazone (5-Chloro-2-Hydroxybenzoxazole) Cholesterol Cimetidine
---	---

Cinchonidine Cinoxacin Ciprofloxacin Citalopram* Citalopram Hydrobromide* Clarithromycin Clemastine Clenbuterol Clindamycin Clindamycin Phosphate Clobazam (except BZO) Clobetasone Butyrate Clonazepam (except BZO) Clonidine Clorazepate(except BZO) Clorazepate Dipotassium Cloxacillin Clozapine Coca ethylene (except COC) Cocaine (Ecgonine Methyl Ester Benzoate) (except COC) Codeine (Desferrioxamine Mesylate) (except BUP, OPI & OXY) Colchicine Cortisone β-Cortol Creatinine Cromolyn (Cromoglycic Acid) Cyclobenzaprine (except TCA) Clomipramine (except TCA) Cyclophosphamide Cyclosporin A Cyproheptadine (except TCA) Dantrolene Deferoxamine Mesylate Deoxyepinephrine R-(-)-Deprenyl (Selegiline) Desipramine (except TCA) N-Desmethylozapine (Normethylozapine) Desmethyldiazepam (except BZO) Desoximetasonone Dexamethasone Dexbrompheniramine Dextromethorphan 4,4'-Diaminophenyl Sulfone (Dapsone) Diazepam (except BZO) Diazoxide Dichloromethane (Methylene Chloride) Dichlorophenamide Diclofenac Dicyclomine Dieldrin Diethyldithiocarbamic Acid N,N-Diethylnicotinamide (Niacin Diethylamide; Nikethamide) Diflorasone Diacetate Diflucortolone pivalate Diflunisal Digitoxin Digoxin (1,2 β-Hydroxydigitoxin) DL-3-4 Dihydroxymandellic Acid DL-3-4 Dihydroxyphenyl Glycol 3,4 Dihydroxyphenylacetic Acid (2,3-Dihydroxypropyl) Theophylline (Dyphylline) Diltiazem Diltiazem-cardzem Dimenhydrinate Dimercaprol (2,3-Dimercaptopropanol) 4-Dimethylaminoantipyrine (Amipopyrine) 1,1-Dimethylbiguanide (Metformin) Dimethyl isosorbide Dimethyl Sulfoxide (DMSO) 1,3-Dimethyluric Acid 1,7-Dimethylxanthine Diphenhydramine (except MTD) 5,5-Diphenylhydantoin (Phenytoin) (except BAR)	Dipyridamole Dipyron Disopyramide (except MTD) Divalproex Dobutamine Doxepin (except TCA) Doxycycline Doxylamine Droperidol Ecgonine Ecgonine Methyl Ester Efavirenz Emetine Enalapril (-)-ψ-Ephedrine (+)-ψ-Ephedrine (+)-Ephedrine (+/-)-Ephedrine (-)-Epinephrine (+/-)-Epinephrine Erythromycin Escitalopram Eserine (Physostigmine) (except OPI) Estazolam (except BZO) β-Estradiol Estril Estrone Estrone-β-D-Glucuronide Estrone-3-Sulfate Ethacrynic Acid Ethambutol Ethamivan (N,N-Diethylvanillamide) Ethanol, Standard Ethopropazine Ethosuximide (2-Ethyl-2-Methylsuccinimide) 2-Ethyl-2-Phenylmalonamide Ethylene Glycol Ethylenediaminetetraacetic Acid (EDTA) 2-Ethylidene-1,5-Dimethyl-3,3-diphenylpyrrolidine Ethylmorphine* (except OPI & OXY) 17-α-Ethinylestradiol Etodolac Etoposide Ezetimibe Famotidine Felodipine Fenfluramine Fenopropfen [(+/-)-2-(3-Phenoxyphenyl) Propionic Acid] Fentanyl* Ferrous Sulfate Fexofenadine Fluoxetine Flurbiprofen Flufenamic Acid Flunisolide Flunitrazepam (except BZO) Fluphenazine Flurandrenolide Flurazepam (except BZO) Flurbiprofen Formaldehyde Furosemide Gabapentin Gemfibrozil Gentamicin Sulfate Gentisic Acid Glucose (D)-(+)-Glucose (Dextrose) Glibenclamide Griseofulvin Guaiacol Glyceryl Ether Guaifenesin Guanethidine Halazepam Halcinonide Haloperidol Hemoglobin	Heroin (Diacetylmorphine)* (except OPI) Hexachlorocyclohexane Hexachlorophene Hexobarbital Hippuric Acid Histamine [2 (4-Imidazolyl) Ethylamine] DL-Homatropine Hydralazine (1-Hydrazinophthalazine) (1S,9R)- β-Hydrastine Hydrochlorothiazide Hydrocodone (except BUP,OPI& OXY) Hydrocortisone Hydroflumethiazide Hydromorphone (except OPI & OXY) Hydroxocobalamin O-Hydroxyhippuric Acid 5-Hydroxyindole-3-Acetic Acid 5-Hydroxy-2-indole-2-Carboxylic Acid 4-Hydroxy-3-Methoxyphenylacetic Acid (Homovanillic Acid) 4-Hydroxy Phencyclidine (except PCP) 11-Hydroxy-Δ9-Tetrahydrocannabinol* (except THC) 5-Hydroxytryptamine (Serotonin) 3-Hydroxytyramine Hydroxyzine (Atarax) L-Hyoscyamine Ibuprofen Irbesartan Imidazole-4-Acetic acid Imipramine (except TCA) Indapamide Indole-3-Acetic acid Indole-3-Butyric Acid DL-Indole-3-Lactic Acid Indomethacin Interferon Ipratropium Bromide Iproniazid Isonicotinic Acid (Pyridine-4-Carboxylic Acid) Isonicotinic Acid Hydrazide Isopropamide (except MTD) (+)-Isoproterenol (-)-Isoproterenol (+/-)-Isoproterenol Isoxsuprine Kanamycin Ketamine Ketoprofen Kynurenin Acid Labetalol Lamotrigine Lanoprazole Lansoprazole Levorphanol Levothyroxine Lidocaine Linoleic Acid-Conjugated (CLA), Gamma, Alpha; Eicosapentaenoic, docahexaenoic acid; omega 369 Lisinopril Lithium Carbonate Loperamide Loratadine (+/-)-Lorazepam (except BZO) Lormetazepam (except BZO) Lysergic Acid Diethylamide (LSD) (except BUP) Mebendazole Meclizine Meclofenamic Acid Medazepam Mefenamic Acid Melanin Meloxicam	Melphalan (-)-Menthhol Meperidine Mephensin Mephentermine Meprobamate 6-Mercaptopurine Mersalyl Acid Mescaline (3,4,5-Trimethoxyphenylethylamine) DL-Metanephine Metaproterenol Metaraminol [(+)-m-Hydroxyphenylpropanolamine] (+/-) Methadone (except MTD) (+) Methamphetamine (Methylamphetamine; d-Desoxyephedrine) (except MDMA & MET) (+/-) Methamphetamine (except MDMA & MET) Methanol, Absolute Methaqualone Methazolamide Methotrimeprazine Methoxamine Methoxamine (S)-6-Methoxy-α-Methyl-2-Naphthalene Acetic Acid (Naproxen) Methoxyphenamine 5-Methoxytryptamine 3-Methoxytyramine 2-Methyl-3-(3,4-dihydroxyphenyl)-DL-Alanine 2-Methyl-3-(3,4-dihydroxyphenyl)-L-Alanine 3,3'-Methylene-bis-(4-Hydroxycoumarin) (Dicumarol) Methylene Blue (+/-) 3,4-Methylenedioxyamphetamine (MDA) (except AMP) (+/-) 3,4-Methylenedioxymethamphetamine (MDMA) (except MET & MDMA) (+/-) 3,4-methylenedioxy-n-ethylamphetamine (MDEA) (except MET & MDMA) 1-Methylhistamine 6 α-Methyl-17 α-Hydroxyprogesterone (Medroxyprogesterone) 6 α-Methylprednisolone (Medrol) Methylphenidate (Ritalin) Methyl Salicylate Methyl Viologen (Gramoxone; Paraquat Dichloride) Meticrane Metoclopramide (except BUP & COC) (+/-)-Metoprolol Metronidazole Mexiletine (except AMP) Mianserin Midazolam Milrinone Minaprine Minocycline Mirtazapine (except BZO) Morphine (except BUP,OPI & OXY) Morphine-3-β-D-Glucuronide (except OPI) Mupirocin Nabumetone Nadolol Nafcillin Nalburphine Nalidixic Acid Nalmefene (except BUP) Nalorphine (except OPI) Naloxone Naltrexone (except BUP) Naphazoline α-Naphthalene Acetic Acid β-Naphthalene Acetic Acid
--	--	---	---

α-Naphthol Neomycin Sulfate Nialamide Nicotinic Acid (Niacin) Nifedipine Nitrazepam (except BZO) Nitrofurantoin Nomifensine 11-Nor-Δ8-Tetrahydrocannabinol-9-Carboxylic Acid* (except THC) 11-Nor-Δ9-Tetrahydrocannabinol-9-Carboxylic Acid* (except THC) 11-Nor-Δ9-THC-9-Carboxylic Acid Glucuronide* (except THC) Norclomipramine (except TCA) Norcocaine Norcodeine (except OPI) Nordoxepin (except TCA) Nordiazepam (except BZO) Norethindrone Norfloxacin DL-Normetanephine Normorphine d-Norpopoxyphene (except PPX) Nortriptyline (except TCA) Noscapine Nylidrin Olmesartan Omeprazole Orotic Acid (Uracil-6-Carboxylic Acid) Orphenadrine Oxalic Acid (Ethanedioc Acid) Oxaprozin Oxazepam (except BZO) Oxolinic Acid Oxybutylin Chloride Oxycodone (except OPI & OXY) Oxymetazoline Oxyphenbutazone Oxyprenolol Oxypurinol Paclitaxel Pancuronium Bromide Pantoprazole Papaverine Pargyline Paroxetine HCL Phenazopyridine Phencyclidine Morpholine (except PCP) Penicillin G (Benzylpenicillin) Pentachlorophenol Pentobarbital (Nembutal) (except BAR) Pentoxifylline (Trental) Pentylene-tetrazole Phencyclidine (except PCP) Phendimetrazine p-Phenylenediamine Phenelzine Phenformin Pheniramine Phenobarbital (except BAR) Phenol Phenolphthalein Phenothiazine (Thiodiphenylamine) Phenoxy-methyl Penicillinic Acid (Penicillin V) Phentermine (α,α-Dimethylphenethylamine) (except AMP) Phentolamine DL-Phenylalanine L-Phenylalanine Phenylbutazone	L-Phenylephrine (+/-)-α-Phenylethylamine (α-Methyl benzylamine) β-Phenylethylamine (R)-(+)-α-Phenylethylamine (+/-) Phenylpropanolamine (PPA) Phenylsamide Phthalic Acid (1,2-Benzenedicarboxylic Acid) Picrotoxin Pilocarpine Pimozide Pinacidil Pindolol Pioglitazone L-Pipecolic Acid Pipemidic Acid Piroxicam Potassium Chloride Potassium Iodide Prazepam Prazosin Prednisolone (1-Dehydrocortisol) Prednisone (Dihydrocortisone) 5-Pregnen-3β-OL-20-one (EPI pregnanolone; Pregnenolone) Prilocaine Primaquine Primidone (2-Desoxyphenobarbital) Proadifen Probenecid [p-(Dipropylsulfamoyl) Benzoic Acid] Procainamide Procaine (Novocaine)(except COC, MDMA & MET) Prochlorperazine Procyclidine (except MTD) Promazine (except TCA) Promethazine Propionyl promazine d-Propoxyphene (except PPX) DL-Propranolol 2-Propylpentanoic Acid (Valproic Acid) Protein Pyridoxine Protriptyline (except TCA) d-Pseudoephedrine Pyridine-2-AldoximeMethochloride (Pralidoxime Chloride) Pyrilamine (Mepyramine) Quinapril Quinidine Quinine Quinolinic Acid (2,3-Pyridinedicarboxylic Acid) Ramipril Ranitidine (Zantac) (except MDMA & MET) Rescinnamine Reserpine Ribavirin Riboflavin Ritodrine Rosiglitazone Rosuvastatin Salbutamol (Albuterol) Salicylamide (2-Hydroxybenzamide) Salicylic Acid (2-Hydroxybenzoic Acid) (-) Scopolamine (Hyoscine) Secobarbital (Quinalbarbitone) (except BAR) Sertraline	Simvastatin Sodium Chloride Sodium Formate (+/-)Sotalol Strychnine Succinylcholine Chloride Sulfamethazine Sulfamethoxazole Sulfanilamide (p-Aminobenzene-sulfonamide) Sulfathiazole Sulfisoxazole Sulindac (except BZO) (+/-)Sulpiride Suxibuzone (except MTD) Talbutal (except BAR) Tamoxifen Tannic Acid Temazepam (except BZO) Tenoxicam Terazonin Terazosin Terazosin HCl Terbutaline Terfenadine Tetracycline Tetraethyl Thiuram Disulfide (Disulfiram) Δ8-Tetrahydrocannabinol (except THC) Δ9-Tetrahydrocannabinol (except THC) Tetrahydrozoline Thebaine (Paramorphine) (except OPI & OXY) Theobromine (3,7-Dimethylxanthine) Theophylline (1,3-Dimethylxanthine) Thiamine (Aneurine) Thimerosal (Sodium Ethylmercurithiosalicylate) Thioridazine cis-Thiothixene Thymol (5-Methyl-2-Isopropylphenol) Timolol Tobramycin Tolazamide Tolbutamide Tolmetin Toluene cis-Tramadol Trans-2-Phenylcyclopropylamine (Tranlycypromine-mine) Tramadol HCl Trazodone Triamcinolone (Fluoxiprednisolone) Triamterene Triazolam* (except BZO) Trichlormethiazide Trichloroacetic acid 2,2,2 Trichloroethanol Trifluoperazine Triflupromazine DL-Trihexyphenidyl Trimethobenzamide (except MDMA & MET) Trimethoprim 3,5,5-Trimethyloxazolidine-2,4-dione (Trimethadione) Trimipramine (except TCA) Triprolidine DL-Tropic Acid Tropine	Tryptamine [3-(2-Aminoethyl) Indole] DL-Tryptophan (3 β-Indolylalanine; (+/-)-α-Amino-3-Indolepropionic Acid) d-Tubocurarine Chloride Tyramine (4-Hydroxyphenethylamine) DL-Tyrosine Urea (Carbamide) Uric Acid Vancomycin (+/-)Verapamil Venlafaxine (except PCP) Vincamine Vitamins Warfarin Xylometazoline Yohimbine Zearalenone Zolpidem Zomepirac Zopiclone
--	--	---	--

* tested at 10,000 ng/mL